



# NEOVELLE<sup>TM</sup> DEFINE



A safe, easy-to-use filler  
ideal for enhancing structure and definition.



Improves facial  
volume!



[ MADE IN ITALY ]

NEOVELLE<sup>TM</sup> DEFINE

## PURE AND MANAGEABLE

**High degree of purity for  
an easy and safe use.**

Neovelle Define is characterized  
by a high degree of purity.

Practical and reliable, its  
formulation is specifically  
designed to deliver a higher  
level of structure and definition.

The patient showed immediate  
improvement with natural results  
and no extended recovery time.



### TECHNICAL DATA SHEET

#### Composition:

- ◆ Cross-linked hyaluronic acid

#### Concentration of hyaluronic acid:

- ◆ 20 mg/ml

#### Packaging:

- ◆ Syringes 2 x 1 ml
- ◆ Needles 4 x 27G 1/2"


## STRENGTHS

- ◆ Ultra-pure hyaluronic acid
- ◆ Enhanced treatment durability and resistance
- ◆ Formulated to restore or improve facial volumes

## APPLICATION

Neovelle Define is injectable into the middle dermis to provide greater structure and definition to the most critical areas of the face. It is also ideal for non-surgical rhinoplasty.

## DERMA INJECTION: MEDIUM



## TREATMENT AREAS

- ① Eyebrow angle lift
- ② Non-surgical rhinoplasty
- ③ Nasolabial fold
- ④ Labial commissure and marionette lines
- ⑤ Periorbital region
- ⑥ Posterior zygomatic region





# CLINICAL CASES



Thanks to the treatment with Neovelle Define, injected into the middle dermis, it is possible to non-surgically correct irregularities and imperfections of the nose, achieving a more harmonious and defined nasal profile.

**BEFORE**



**AFTER**



Definition of the posterior zygomatic region



**BEFORE**



**AFTER**



Non-surgical rhinoplasty



**ESTIMATED DURATION:  
UP TO 6 MONTHS**

Treatment: Zygomatic area  
Amount: 1 ml  
Filling capacity: ●●○○

## PATIENT RESULTS:

**98%**

Of patients report an improvement in their appearance

**95%**

Find their results to be natural

**BEFORE**



**AFTER**



Correction of the labial commissure

**BEFORE**



**AFTER**



Definition of the posterior zygomatic region and the nasolabial fold

# NEOVELLE™

Extra Pure Hyaluronic Acid



**High-quality, cross-linked and highly purified hyaluronic acid. Designed to restore skin volume while preserving natural contours.**

Neovelle is a range of dermal fillers formulated with the purest hyaluronic acid, available in 3 specific variants, for targeted results across different areas of the face. Designed to restore skin volume while preserving natural contours.

Our fillers contain a BDDE (butanediol diglycidyl ether) residue of less than 0.08 ppm, do not cause edema, have a high corrective capacity, and provide a smooth, even, and full result.





“



**Marllor Biomedical has been a global leader in creating and distributing medical devices and procedures for aesthetic medicine, combining innovation, uniqueness, and attention to detail.**

”

Marllor Biomedical's mission is to put our 15 years of experience at the service of aesthetic medicine and dermatology, offering innovative solutions tailored to specialists in the field. We are dedicated to providing healthcare professionals with effective, safe and targeted treatment solutions that address patients' needs. In the field of non-surgical aesthetic medicine, Marllor Biomedical is synonymous with innovation, quality and attention to detail on a global scale.



MADE IN ITALY



## ✧ NEOVELLE<sup>TM</sup> SOFT

The most delicate and natural filler in the Neovelle line, it reduces fine lines and enhances lip volume, providing a hydrating effect that lasts up to 6 months.



## ✧ NEOVELLE<sup>TM</sup> DEFINE

Safe and precise structural correction, perfect for addressing deep wrinkles and redefining facial contours. The high purity and viscosity ensure harmonious and long-lasting results.



## ✧ NEOVELLE<sup>TM</sup> VOLUME

The densest and most volumizing filler in the range, perfect for areas that require added structure and support. It provides visible, natural-looking results that last up to 10 months.

INFORMATION RESERVED FOR THE MEDICAL CORPORATE. DISCLOSURE TO THE PUBLIC IN ANY FORM IS PROHIBITED

### MARLLOR BIOMEDICAL

#### REGISTERED OFFICE

Via Torconca, 1 - 47842  
S. Giovanni in Marignano (RN)  
Partita IVA 03189570405

☎ +39 0541 410286  
✉ [info@marllor.com](mailto:info@marllor.com)  
🌐 [www.marllor.com](http://www.marllor.com)

SCAN QR CODE

